Miltenyi Biotec

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REAfinity™ Antibodies

ISOTYPE CONTROL

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Realease Fluorochrome Technology

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
### Manufacturing workflow for REAfinity Antibodies

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



Mass spectrometry analysis of purified antibodies shows that REAfinity Antibodies are defined products, while hybridoma generated antibodies can be a mixture. (A) Mass spectrometry analysis of two examples of hybridoma-generated monoclonal antibodies shows that both contain a second light chain with a molecular weight of approx. 23,635 Da. In one example bottom left, the amount of contaminating light chain was by far exceeding the productive light chain of approx. 24,142 Da. (B) Two examples of recombinant REAfinity Antibodies show pure light chain populations.

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Human PBMCs from a single donor were stained with CD56 antibodies conjugated to FITC, PE, APC, APC-Vio® 770, or PE-Vio® 770 from two different production lots, and analyzed on the MACSQuant® Analyzer. The overlayed histograms represent the different lots. Tandem Signal Enhancer, human, was used to increase stain indices of tandem dye-conjugated antibodies. Cell debris and dead cells were excluded from the analysis based on scatter signals and PI fluorescence. In case of staining with CD56 conjugated to tandem dyes, DAPI was used to exclude dead cells from the analysis.

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
Other clones	Overlap in epitope recognition with REA196
NCAM16.2	++
TULY56	++
5.1H11	++
AF12-7H3	-
C5.9	++
N901/NKH1	++
CMSSB	+
B159	-
HCD56	-
MEM188	-

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
Fluorescence microscopy image of CD53 knockout cells. Wild type (WT, left) and knockout cells (KO; right) were stained with CD53-PE (REA259, red) and counterstained with DRAQ5 (blue) as DNA stain.

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Knockdown of the target antigen by RNAi. The left peak indicates HeLa cells transfected with AKT1 RNAi, the right peak is the non-transfected control. Cells were fixed, permeabilized, stained with Anti-AKT1 antibody-fluorochrome conjugates, and analyzed by flow cytometry using the MACSQuant Analyzer 10. Cell debris was excluded based on scatter signals.

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

High staining index with REAfinity Recombinant Antibodies. Human PBMCs were stained with different hybridoma clones and REAfinity Clone REA196 against CD56. Recommended titer of respective antibodies from different suppliers was used.

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Comparison data for clone REA196. Performance of fluorochrome-conjugated CD56 antibodies (clone REA196) were compared to two commercially available clones from competitors in human PBMCs on a MACSQuant Analyzer 10.

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
\* Exclusion channel containing lineage markers to exclude T cells (CD3) and B cells (CD19) from the analysis.

#### Gating strategy

Splenocytes from C57BL/6 mice were stained using the described antibody panel to analyze NK cells. Mouse NK cells are identified within the CD45+/lineage (CD3/CD19) negative population as NK1.1+ and CD335+ cells. Expression of KLRG1 on mouse NK cells correlates with mature NK cells. NK cells can be divided into four maturation stages based upon the expression of CD27 and CD11b. Stages in order from less to more mature are: S1 < S2 < S3 < S4.

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

Splenocytes from BALB/c mice were stained with a PE-conjugated Anti-MHC Class II antibody (B) or with the corresponding Rat IgG2b antibody (A).

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### FcR blocking reagents

130-059-901	FcR Blocking Reagent, human	2 mL
130-092-575	FcR Blocking Reagent, mouse	2 mL

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
Human whole blood samples (1 mL) from a SARS-CoV-2-reactive donor were incubated for 8 hours with Brefeldin A and the SARS-CoV-2 PepTivator Prot\_N (130-126-700) or Prot\_S1 (130-127-041), or left unstimulated (negative control). Blood samples were lysed, fixed, and permeabilized. Subsequently, cells were stained with the flow panel included in the SARS-CoV-2 T Cell Analysis Kit (Whole Blood), human. Cells were analyzed using a MACSQuant Analyzer 16. Doublets, debris, and dead cells as well as CD14+ and CD20+ cells were excluded. After pre-gating on CD3 as well as CD4 and CD8, respectively, activation marker and cytokine expression were assessed, e.g., CD154 and IL-2 for CD4+ T cells and TNF- $\alpha$  and IFN- $\gamma$  for CD8+ T cells.

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Marker	Expressed by
MSC-positive marker	
CD73, CD90, CD105	MSCs
MSC-negative marker	
CD14	Monocytes and macrophages
CD19	Pan B cells
CD34	Primitive hematopoietic progenitors and endothelial cells
CD45	Pan-leukocyte cells
HLA-DR	Dendritic cells, B cells, monocytes, macrophages

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
### Myeloid-derived suppressor cells (MDSCs) 細胞標記

M-MDSC	HLA-DR-CD11b+CD14+CD33 <sup>high</sup>
PMN-MDSC	HLA-DR-CD11b+CD14-CD33 <sup>Mid</sup>

### MDSC Detection Kit human

130-126-233	50 tests
1 mL MDSC Staining Cocktail, human 包含以下抗體	
7-AAD /CD11b-VioBlue (clone: REA713)/CD45-VioGreen™ (clone: REA747)/HLA-DR FITC(clone: REA805)/CD33-PE (clone: REA775)/CD16-PE-Vio® 770(clone: REA423) /CD14 -APC-Vio 770(clone: REA599)	
1 mL MDSC Control Cocktail, human 包含以下抗體	
7-AAD and CD11b -VioBlue (clone: REA713)/ CD45-VioGreen (clone: REA747)/ REA Control (S)-FITC(clone: REA293)/REA Control (S)-PE(clone: REA293)/CD16-PE-Vio 770(clone: REA423)/CD14-APC-Vio 770(clone: REA599).	
預留 APC 偵測通道	
可依據實驗需求加入 CD66b-APC, human (clone: REA306), CD15-APC, human (clone: VIMC6), Anti-LOX-1-APC, human (clone: REA1188), CD10-APC, human (clone: REA877), REA Control (S)-APC.	
0.1 mL CD16 antibody, anti-human, PE-Vio 770, REAfinity:for compensation control	
0.1 mL CD14 antibody, anti-human, APCVio 770, REAfinity:for compensation control	

### Mouse MDSC detection maker

### Mouse

品項	clone
Ly-6C Antibody, anti-mouse, REAfinity™	REA796
CD115 Antibody, anti-mouse, REAfinity™	REA827
CD192 (CCR2) Antibody, anti-mouse, REAfinity™	REA538
Ly-6G Antibody, anti-mouse, REAfinity™	REA526
CD11b Antibody, anti-mouse, REAfinity™	REA592

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RUO=Research Use Only

\* For complete regulatory and legal notices please refer to the back page of the brochure.

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## REAfinity™ Antibodies for mass cytometry

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

#### Comparison between flow and mass cytometry.

(A) Human peripheral blood mononuclear cells (PBMCs) were stained with CD159c (NKG2C)-PE (REA205) and analyzed using flow cytometry

(B) The same sample was stained with the metal-conjugated pure version of CD159c (NKG2C) (REA205) and analyzed using mass cytometry. Data courtesy of Prof. Karl Johan Malmberg, Oslo University Hospital, Norway.

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#### Expression pattern of NKG2A and NKG2C on NK cells.

(A) Representative mass cytometry plots show NK cell marker expression including NKG2A and NKG2C counterstained with CD56.


(B) Additionally, viSNE analysis was performed by clustering 20,000 NK cells according to the expression of 13 different markers. The following REAfinity Recombinant Antibodies were used for staining: CD56 (REA196), CD122 (IL-2R  $\beta$ ) (REA167), CD38 (REA671), CD39 (REA739), CD366 (TIM-3) (REA635), TIGIT (REA1004), CD226 (DNAM-1) (REA1040), CD159a (NKG2A) (REA110), CD159c (NKG2C) (REA205), CD158a (KIR2DL1) (REA284), CD158a/h (KIR2DL1/DS1) (REA1010), CD158b (KIR2DL2/DL3) (REA1006), and CD158e (KIR3DL1) (REA1005). Data courtesy of Dr. Amir Horowitz, Icahn School of Medicine at Mount Sinai, United States.

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# REAl ease® Fluorochrome Technology

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
## Automated features of the MACSQuant Family

Startup, cleaning,  
and shutdown

Compensation

Labeling

Calibration and  
quality control (QC)

## Express Mode Master Package (Cat. No. 160-002-425)

CD19 CAR T Cell Express Mode Package	
Immune_Cell_Composition_h_02	鑑定免疫細胞存活率與組成
CAR_T_Cell_Transduction_h_02	鑑定 CART 轉導效率
B_Malignant_Cells_h_01	鑑定惡性 T 細胞
B_Regulatory_Cells_h_01	鑑定調節型 B 細胞亞群
B_Cell_Differentiation_h_02	鑑定 B 細胞分化
CAR_T_Cell_Tscm_h_01	鑑定 stem memory T cells (Tscm)
CAR_T_Cell_Persistence_h_01	鑑定 CART 持續性
CAR_T_Cell_Differentiation_h_01	鑑定 CAR+ and CAR- T cells 細胞分化
CAR_T_Cell_Exhaustion_h_01	鑑定 CAR+ and CAR- T cells 細胞衰老
CAR_T_Cell_Activation_h_01	鑑定 CAR+ and CAR- T cells 活化
CAR_T_Cell_Staining_Control_h_01	染色對照組
TCR $\alpha$ $\beta$ /CD45RA Express Mode Package	
TCR $\alpha$ $\beta$ _CD45R_depletion_h_02	鑑定 CD45RA+, TCR $\alpha$ / $\beta$ +, and CD34+ cells 比例
Immune_Cell_Composition_human	鑑定免疫細胞存活率與組成
Virus-Specific T Cell CCS Express Mode Package	
RCI_CD4CD8_h_02	鑑定 IFN- $\gamma$ 表達的目標細胞
CCS_Immune_Cell_Composition_h_01	鑑定細胞組成與目標細胞比例
CCS_Purity_h_01	鑑定目標細胞純度

## CAR T Cell Essential Express Mode Package (Cat. No. 160-002-759)


Express mode	Purpose	V1 VioBlue®	V2 VioGreen™	B1 FITC	B2 PE	B3 7AAD Per-CP- Vio® 700	B4 PE- Vio® 770	R1 APC	R2 APC- Vio® 770
CAR T cell immune composition	For the determination of cell viability and composition of immune cells	CD45	CD4	CD3	CD56/ CD16	7AAD	CD19	CD14	CD8
CAR T cell transduction	For the determination of transduction efficiency of CAR T cells	CD45	CD4	CD3	CAR detection reagent	7AAD	-	CD14	CD8
CAR T cell persistence	For the determination of CAR T cell persistence	CD45RA	CD4	CD3	CAR detection reagent	7AAD*	CD62L	CD45RO	CD8

\*For immunomonitoring, exclusion markers such as CD14 and CD15 conjugated to PerCP-Vio® 700 can be added

## CAR T Cell Essential Express Mode Package

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Cat.no.	Name	Size
130-115-965	CD19 CAR Detection Reagent, human, Biotin	10 tests in 100 µL
130-126-727	CD22 CAR Detection Reagent, human, Biotin	30 tests in 60 µL
130-126-090	BCMA CAR Detection Reagent, human, Biotin	30 tests in 60 µL
130-113-853	Biotin Antibody, PE	30 tests in 60 µL
130-113-291	Biotin Antibody, PE	100 tests in 200 µL

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