


KAPA3G Plant PCR Kits




做植物 PCR
終於有更好的選擇啦！

- 簡便省時**：不需抽取植物 DNA，直接進行 PCR
- 高成功率**：對植物抑制劑 (如：**polyphenolics**) 有高耐受性，PCR 成功率高
- 速度快**：可應用於大量轉殖基因 (**GMO**) 植物的檢測
- 靈敏度高**：只需非常少量的組織，且不需研磨成粉末
- 應用性廣**：PCR 產物還可以做 **T/A Cloning**



Cat No	Product	Size
KK7251	KAPA3G Plant PCR Kit	250 x 50 μ L rxn
KK7252	KAPA3G Plant PCR Kit	500 x 50 μ L rxn


從葉片、種子和其他植物組織中直接進行 PCR



各種植物、不同組織皆可獲得良好結果

PCRs (50 μ L) contained the crude sample or 1-10 ng purified DNA (depending on the species), and 40 cycles were performed in all cases. Targets ranged between 500 and 900 bp, and reaction products were analyzed in a 1% agarose gel. KAPA Express DNA Ladder (100, 200, 400, 800, 1600, 4000, 8000 bp) was used as a MW marker.


大量節省您寶貴的實驗時間



傳統方法 (CTAB) 耗時、費力 (萃取 DNA+PCR 約需 3.5 小時)，KAPA3G Plant PCR Kit 僅需 45 分鐘即可完成 PCR。

CTAB-extracted DNA and crude samples were used as templates in 50 μ L PCRs, with 40 cycles of amplification. Direct PCR with the KAPA3G Plant PCR Kit was completed in 45 min. In contrast, the CTAB extraction protocol required ~2 h, and the amplification with wild-type Taq 1.5 h to complete.

KAPA3G Plant PCR Kit 較其他廠牌 Taq 具有更高靈敏度



KAPA3G Plant PCR Kit 與它牌的試劑結果比較

不論純化過的 DNA 或是粗粹取的樣本成功率均勝出

Purified DNA (1-10 ng per reaction, depending on species) and crude samples were used as templates in 50 μ L PCRs, with 40 cycles of amplification. Targets ranged between 500 bp and 900 bp and reaction products were analyzed in a 1% agarose gel.