Miltenyi Biotec

# TumorMACS™ Media

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### Collection and shipment of samples

- Sample collection
- Storage of tissue in MACS® Tissue Storage Solution
- Shipment, if necessary


### Tumor dissociation

- Tumor Dissociation Kit, human for primary human or human xenografted tumors
- Cultivation\* of unsorted cells (e.g. including fibroblasts)

### Isolation of tumor cells

- Tumor Cell isolation Kit, human or Mouse Cell Depletion Kit for isolation of tumor cells
- Cultivation\* of pure tumor cells

### Flow cytometric analysis

- Determination of sample composition
- Phenotyping of tumor cells

### \*Cultivation


- Plate  $2 \times 10^5$  cells/cm<sup>2</sup> in TumorMACS™ Medium
- Incubate at 37 °C and 5–7.5% CO<sub>2</sub>
- Change medium every 4–5 days or as soon as it turns yellowish
- Split or remove fibroblasts at confluency of 80–90%

Workflow overview on how to establish a stable cell line from primary tumor cells.

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
PANC-1-derived xenografts show a homogeneous histology.

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
Patient-derived pancreatic tumor cells were injected into mice. Tumor tissue was extracted from the xenograft, cultured in Pancreas TumorMACS Medium over multiple passages, and re-injected into mice for the validation of tumorigenicity.

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### Establishment of adherent or suspension culture.

The images show two different cell lines, derived from the same primary pancreatic adenocarcinoma, growing in adherent or suspension culture. During initial cultivation of the primary tumor cells, some cells detached, grew in suspension (B), and gave rise to a separate cell line, distinct from the cells that kept their adherent characteristics (A). Pancreas TumorMACS Medium supported cell expansion both in adherent and suspension cultures. Remarkably, in the samples shown, the epithelial tumor cells grew better in a suspension culture of spheroid-like structures (B vs. A).