# NEXT GENERATION SEQUENCING

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



# KAPA HYPERPLUS GUIDE TO SUCCESS

## Got EDTA?

- The enzymatic fragmentation reaction is sensitive to EDTA.
- The best strategy is to remove EDTA by means of a cleanup step before fragmentation.
- If your DNA contains EDTA, please see Appendix 2 (p. 16) of the Technical Data Sheet (TDS).

## How much DNA do I need?

Application	Sample type	Recommended input
WGS	Complex gDNA (high quality)	50 ng – 1 µg
Target capture (WES, custom panels)	Complex gDNA (high quality)	10 ng – 1 µg
WGS, target capture	FFPE DNA	≥50 ng (quality dependent)
WGS	Microbial DNA	1 ng – 1 µg
WGS (PCR-free)	High-quality DNA	≥50 ng (no SS)* ≥500 ng (w/SS)*
Targeted sequencing	Long amplicons	≥1 ng

\*SS = size selection

## Get to chopping.

- Mode and size distribution of DNA is controlled by fragmentation time and temperature.
- Try a range of fragmentation times to determine optimal insert size.
- For ease of sample processing, place samples with the longest fragmentation time in the thermal cycler first. Add samples with shorter fragmentation times at appropriate intervals.

Mode fragment length	Incubation time at 37°C*	Optimization range
600 bp	5 min	3 – 10 min
350 bp	10 min	5 – 20 min
200 bp	20 min	10 – 25 min
150 bp	30 min	20 – 40 min

\*These parameters are a good starting point for high-quality genomic DNA. Please refer to Appendix 2: Optimization of Fragmentation Parameters of the TDS for guidelines on how to optimize fragmentation time and temperature, if needed.

## It's not a typo!

- Ensure that you are adding the correct volume of KAPA Frag Buffer (5 µL) and KAPA Frag Enzyme (10 µL) to each reaction.

Component	Volume
Double-stranded DNA (with KAPA Frag Conditioning Solution, if needed)	35 µL
KAPA Frag Buffer (10X)	5 µL
KAPA Frag Enzyme	10 µL
<b>Total volume</b>	<b>50 µL</b>

## How much adapter

- Adapter concentration as adapter and adapter-ligation cleanup.

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
KAPA HyperPlus: single-tube fragmentation and library preparation workflow in less than 3 hours

### 4 ADAPTER LIGATION

incubate at 20°C

15 min

**How much do I need?**  
 affects ligation efficiency, as well as adapter dimer carry-over during the post-ligation cleanup.

Adapter stock concentration	Adapter:insert molar ratio
15 μM	10:1
15 μM	20:1
15 μM	40:1
15 μM	100:1
15 μM	200:1
7.5 μM	200:1
3 μM	200:1
1.5 μM	200:1
750 nM	200:1
300 nM	200:1



### Don't do more work than you have to.

- Determine how much final library material is required for your downstream application (e.g., capture, sequencing).

Input into library construction	Number of cycles required to generate	
	100 ng library	1 μg library
1 μg	0	0 - 1
500 ng	0	2 - 3
250 ng	0 - 1	3 - 5
100 ng	0 - 2	5 - 6
50 ng	3 - 5	7 - 8
25 ng	5 - 6	8 - 10
10 ng	7 - 9	11 - 13
5 ng	9 - 11	13 - 14
2.5 ng	11 - 13	14 - 16
1 ng	13 - 15	17 - 19

### So did you make a good library?

This is where you decide if you made a good library by assessing both the size and concentration with an appropriate quantification method.

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

The ligation of adaptor sequences (orange and blue) to library DNA molecules (grey) results in a mixed population of molecules without the correct adaptor configuration. Electrophoresis and spectrophotometry measure total nucleic acid concentrations, whereas optimal cluster density or template-to-bead ratio depends on the appropriate concentration of PCR-amplifiable DNA molecules.

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
The KAPA Library Quantification Kit was used to determine the concentration of two Illumina GA libraries with unusual GC content (*Rhodococcus* sp.; ~70% GC - shown above, *Staphylococcus* sp.; ~35% GC - not shown). Both libraries amplified with efficiency >95%. Two-fold dilution series (1:1000 through 1:16000) were prepared in triplicate, and qPCR performed according to the recommendations in the product technical data sheet.

Cat No.	Product	Size
07960140001	Complete kit (Universal)	500 x 20 µl rxns
07960204001	Complete kit (ABI Prism®)	500 x 20 µl rxns
07960255001	Complete kit (Bio-Rad®)	500 x 20 µl rxns
07960298001	Complete kit (optimized for Roche® LightCycler 480)	500 x 20 µl rxns
07960336001	Complete Kit (ROX Low)	500 x 20 µl rxns

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
Fragmented DNA per 50 $\mu$ L ER & AT reaction	Adapter stock concentration	Adapter:insert molar ratio	Fragmented DNA per 50 $\mu$ L ER & AT reaction	Adapter stock concentration	Adapter:insert molar ratio
1 $\mu$ g	15 $\mu$ M	10:1	25 ng	7.5 $\mu$ M	200:1
500 ng	15 $\mu$ M	20:1	10 ng	3 $\mu$ M	200:1
250 ng	15 $\mu$ M	40:1	5 ng	1.5 $\mu$ M	200:1
100 ng	15 $\mu$ M	100:1	2.5 ng	750 nM	200:1
50 ng	15 $\mu$ M	200:1	1 ng	300 nM	200:1

\* Adapter:insert molar ratio calculations are based on a mode DNA fragment length of 200 bp, and will be higher for longer DNA fragments, or slightly lower for DNA fragmented on a mode size <200 bp. The lower adapter:insert molar ratios recommended for inputs >100 ng represent a fair compromise between library construction efficiency and cost; higher library yields will be achieved if a higher adapter concentration is used.

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

Libraries were prepared from 50 ng hgDNA or 50 ng FFPE DNA and captured with the Nimblegen™ SeqCap EZ HGSC VCRome panel.

Cat No	Product	Size
07962380001	KAPA HyperPlus Kit	8 rxns
07962401001		24 rxns
07962428001		96 rxns
08005702001	KAPA Single-Indexed Adapter Set A (30 $\mu$ M)	12 x 40 $\mu$ L
08005729001	KAPA Single-Indexed Adapter Set B (30 $\mu$ M)	12 x 40 $\mu$ L

for research use only, not for clinical purpose

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



Libraries were prepared from 100 ng of FFPE DNA, and exome capture was performed with the SeqCap EZ HGSC VCRome Design, Roche NimbleGen.

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Library diversity and duplication rates for libraries prepared from 2 ng cell-free DNA, using the KAPA Hyper Prep Kit or a leading competitor kit optimized for low-input library construction.

Cat No	Product	Size
07962312001	KAPA Hyper Prep Kit for Illumina	8 rxns
07962347001		24 rxns
07962363001		96 rxns
08005702001	KAPA Single-Indexed Adapter Set A (30 μM)	12 x 40 μL
08005729001	KAPA Single-Indexed Adapter Set B (30 μM)	12 x 40 μL

for research use only, not for clinical purpose