

全方位純化試劑系列

DNA Clean & Concentrator

- 最快可以在兩分鐘內純化並濃縮 DNA
- 回溶體積最小可達 6 μ l，濃縮純化高品質的 DNA
- 純化範圍涵蓋大部分常用 DNA 片段大小
- 特殊 column 設計，完美回收樣品不殘留

Technology Overview: DNA Clean & Concentrator®

Select-a-Size DNA Clean & Concentrator®

High-quality, size selected DNA in 7 minutes (library preparation and NGS applications)

Genomic DNA Clean & Concentrator® Kits

High molecular weight DNA clean-up (1 kb to > 200 kb)

DNA Clean & Concentrator® Kits


Ultra pure DNA in 2 minutes (50 bp to 23 kb) from PCR, impure preps and enzymatic digestions

Oligo Clean & Concentrator™ Kits

DNA & RNA oligos and probes (16 to 200 nt)

DCC 產品 selection guide

依核酸長度不同選擇產品



Cat. No.	Product	Size	Binding Capacity	Minimum Elution
D4013 系列, D4033 系列, D4029 系列, D4031 系列	DNA Clean & Concentrator-5/25/100/500	50 preps	5/25/100/500 μ g	$\geq$ 6/25/150 μ l/2 mL
		200 preps		
D4010, D4011 D4064, D4065	Genomic DNA Clean & Concentrator-10/25	25 preps	10/25 μ g	$\geq$ 10/35 μ l
		100 preps		
D4060, D4061	Oligo Clean & Concentrator	50 preps	10 μ g	$\geq$ 6 μ l
		200 preps		
D4080	Select-a-Size DNA Clean & Concentrator	25 preps	3 μ g	$\geq$ 10 μ l
D4007, D4008	Zymoclean Gel DNA Recovery Kit	50 preps	5 μ g	$\geq$ 6 μ l
		200 preps		